



# SAFETY DATA SHEET

US Department of Labor  
Occupational  
Safety and Health Administration

Product Name: AQUASONIC CLEAR® ULTRASOUND  
GEL

## SECTION 1 — CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

**Product Name:** AQUASONIC CLEAR® ULTRASOUND GEL  
**(As used on label and List)**

**Manufacturer's Name:** PARKER LABORATORIES INC.  
**Address:** 286 ELDRIDGE ROAD  
**(number, street, City, State, Zip code)** FAIRFIELD, NJ 07004 USA  
**Emergency Telephone Number:** (800)-631-8888  
**Telephone Number for Information:** (973)-276-9500 (8:30am – 5:15pm ET)

## SECTION 2 — HAZARDS IDENTIFICATION

**Hazardous Components (Specific Chemical Identity; Common Name(s)):** NONE  
**Classification of the Substance/Mixture:** NOT CLASSIFIED  
**Signal Word:** N/A  
**Hazard Statements:** NO KNOWN SIGNIFICANT EFFECTS OR CRITICAL HAZARDS  
**Precautionary Statement:** HANDLE IN ACCORDANCE WITH GOOD INDUSTRIAL HYGIENE PRACTICES  
**OSHA Hazards:** NO OSHA HAZARDS

## SECTION 3 — COMPOSITION/INFORMATION ON INGREDIENTS

Hazardous Ingredients (specific)	% Composition	CAS Number
N/A	N/A	N/A

## SECTION 4 — FIRST AID MEASURES

**Signs and Symptoms of Exposure:** N/A  
**Medical Symptoms Generally Aggravated by Exposure:** NONE KNOWN  
**Emergency and First Aid Procedures:**

**Inhalation:** Move to fresh air immediately. If experiencing difficulty breathing, seek medical attention

**Skin contact:** In the case of skin irritation or allergic reactions, wash with soap and plenty of water. See a physician.

**Eye contact:** In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

**Ingestion:** If accidentally swallowed obtain immediate medical attention. Clean mouth with water and drink afterwards plenty of water. Do not induce vomiting without medical advice. Never give anything by mouth to an unconscious person.



# SAFETY DATA SHEET

Product Name: AQUASONIC CLEAR® ULTRASOUND GEL

US Department of Labor  
Occupational  
Safety and Health Administration

## SECTION 5 — FIRE FIGHTING MEASURES

Flash Point(Method Used): NON-FLAMABLE      Flammable Limits: N/A      LEL: N/A      UEL: N/A

Extinguishing Media: USE EXTINGUISHING MEDIA APPROPRIATE FOR SORROUNDING FIRE

Special Fire Fighting Procedures: NONE

Unusual Fire and Explosion Hazards: NONE

Protective Equipment and Precautions for Firefighters: WEAR A SELF-CONTAINED BREATHING APPARATUS AND FULL PROTECTIVE GEAR.

## SECTION 6— ACCIDENTAL RELEASE MEASURES

Leak and Spill Procedures: SLIPPERY IF SPILLED ON FLOOR. MIX WITH SWEEPING COMPOUND. MOP AREA WITH WATER.

Personal Precautions: N/A

Environmental Precautions: N/A

## SECTION 7— HANDLING AND STORAGE

Handling procedures and Equipment: N/A

Storage requirements: NO SPECIAL PRECAUTIONS REQUIRED.

## SECTION 8— EXPOSURE CONTROL/PERSONAL PROTECTION

Exposure Limits:     ACGIH TLV                       OSHA PEL                       Other(Specify)

Specific Engineering Controls (such as ventilation, enclosed process)      ENSURE ADEQAUTE VENTILATION

Personal Protective Equipment:

Gloves    Eye     Clothing     Footwear    Respirator    Other (specify)

If checked, please specify type:

## SECTION 9— PHYSICAL AND CHEMICAL PROPERTIES

Physical State(Appearance): Viscous, Clear Aqueous Gel

Color: Colorless

Odor: None

Specific Gravity(H<sub>2</sub>O=1): 1.02

pH: 6.50-6.95

Solubility in Water: Soluble in water

Boiling point: N/D

Melting Point: N/D

Flammability(solid, gas): N/D



# SAFETY DATA SHEET

US Department of Labor  
Occupational  
Safety and Health Administration

Product Name: AQUASONIC CLEAR® ULTRASOUND  
GEL

Vapor Pressure(mm Hg): N/D  
Vapor Density(AIR=1): N/D  
Viscosity: 130,000-185,000cps  
Density: N/D  
Evaporation Rate(Butyl  
Acetate= 1): N/D

## SECTION 10 — STABILITY AND REACTIVITY

Stability: Unstable       Conditions to Avoid: NONE  
Stable

Hazardous Decomposition or Byproducts: NONE

Incompatibility (Materials to Avoid): NONE

Hazardous Polymerization:      May Occur?       Will Not Occur?

Conditions to Avoid: NONE

Reactivity: NO SPECIFIC TEST DATA RELATED TO REACTIVITY AVAILABLE FOR THIS PRODUCT.

## SECTION 11 — TOXICOLOGICAL INFORMATION

Health Hazards (Acute and Chronic): NONE KNOWN

Route(s) of Entry:  Skin Contact     Eye Contact     Inhalation       Ingestion

Carcinogenicity:      NTP? NO      IARC Monographs? NO      OSHA Regulated: NO

Signs and Symptoms of Exposure: N/A

Medical Conditions Generally Aggravated by Exposure: NONE KNOWN

## SECTION 12 — ECOLOGICAL INFORMATION

Ecotoxicity: N/A  
Persistence and Degradability      THERE IS NO DATA AVAILABLE FOR THIS PRODUCT  
Bioaccumulative Potential      THERE IS NO DATA AVAILABLE FOR THIS PRODUCT  
Mobility in the Soil:      THERE IS NO DATA AVAILABLE FOR THIS PRODUCT  
Additional Ecological Information:      THERE IS NO DATA AVAILABLE FOR THIS PRODUCT

## SECTION 13 — DISPOSAL CONSIDERATIONS

Waste Disposal:

FOLLOW ALL FEDERAL, STATE AND LOCAL REGULATIONS FOR NON-HAZARDOUS WASTE  
DISPOSAL.



# SAFETY DATA SHEET

Product Name: AQUASONIC CLEAR® ULTRASOUND GEL

US Department of Labor  
Occupational  
Safety and Health Administration

## SECTION 14 — TRANSPORT INFORMATION

**DOT Classification:** NOT REGULATED  
**IMDG Classification:** NOT REGULATED  
**TDG Classification:** NOT REGULATED  
**IATA Classification:** NOT REGULATED  
**RID:** NOT REGULATED

## SECTION 15 — REGULATORY INFORMATION

**OSHA Hazards:** NO OSHA HAZARDS  
**SARA 311/312 Hazards:** NO SARA HAZARDS

## SECTION 16 — OTHER INFORMATION

**Date of Issue:** 5/27/2015  
**Prepared By:** Quality Assurance  
**Version:** 1

***Disclaimer:** This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. No warranty, guarantee or representation is made as to the correctness and sufficiency of the information. It's the user's responsibility to determine safe conditions for the use of this product and to assume liability for loss, injury, or expense resulting from the misuse of this product.*