

SANI-CLOTH ACTIVE

Canister	SANI – CLOTH ACTIVE	XP00190	200 COUNT	6 x 200 wipes per case
Bucket	SANI – CLOTH ACTIVE	XP00193	225 COUNT	1 X 225 SU
Canister	SANI – CLOTH ACTIVE Canister	PDI Code XP00188	125 COUNT	6 x 125 wipes per case



Description : Sani-Cloth Active is an alcohol Free disinfectant wipe with proven effectiveness against over 70 micro-organisms. 99.99% effective against bacteria, Fungi, Viruses and TB and Polio

Wipe Size and Pack Formats :

- XP00190 – 200 count 200mmx 220mm – 23gsm material
- XP00188 - 125 count – 133mm x 220mm – 23gsm material
- XP00193 – 225 count – 300mm x 245mm – 25gsm material

Benefits :

- Alcohol Free so less likely to cause damage to surfaces
- Large or small wipe sizes available
- Large canister with Deep Well Lid to ease dispensing
- Flow Wrap with Solo Lid to ease dispensing
- Already used by GE and other companies for Ultrasound probes

Cost Savings

- Small canister wipes for small areas like Ultrasound probes
- Economical Flow Wrap format available

Infection Control:

- 99.99% bacteria, fungi and virucidal kill
- Meets EN1040, EN1276, EN13697, EN1275 & EN14476
- Even effective against non-envelope viruses
- Kills TB, Polio and Bird Flu